


# OPERATING INSTRUCTIONS

## AND SUGGESTED ACTIVITIES

### EQUAL VOLUME METAL CYLINDERS SET



#### DESCRIPTION

The Equal Volume Metal Cylinders Set consists of four identically-sized cylinders of different metals and a storage box with a fitted compartment for the samples. Each cylinder is approximately 50 mm (2") long and 12.7 mm (1/2") in diameter. The metals are aluminum, brass, copper, and steel.

#### BACKGROUND

The cylinders are used for measuring and comparing the densities of the four different metals. They can be used in three different ways:

- If students are not expected to calculate the volume of a cylinder, the four specimens can simply be weighed individually. Since the cylinders' volumes are all equal, the ratio of their weights gives the ratio of their densities.
- The weights and dimensions of the cylinders can be measured, the volume calculated, and the densities determined as "weight per unit volume."
- Using a balance with a specific gravity capability, each cylinder can be weighed in air and in water, and the specific gravity calculated from Archimedes' Principle.

#### SUGGESTED ACTIVITIES

##### ACTIVITY 1

##### COMPARISON OF RELATIVE DENSITIES

##### EQUIPMENT NEEDED

1	Equal Volume Metal Cylinder Set
1	Mechanical or electronic balance capable of measuring to 0.1 g

PROCEDURE

- Set up the balance on a flat, stable surface.
- Weigh each of the cylinders in turn and enter their weights into the table below:

Metal	Weight (grams)
Aluminum	
Brass	
Copper	
Steel	

Which metal sample weighed the most? \_\_\_\_\_

Which metal sample weighed the least? \_\_\_\_\_

Which metal do you think is the most dense? \_\_\_\_\_

Which is the least dense? \_\_\_\_\_

ACTIVITY 2

CALCULATING THE DENSITIES

EQUIPMENT NEEDED

- 1 Equal Volume Metal Cylinder Set
- 1 Mechanical or electronic balance capable of measuring to 0.1 g
- 1 Vernier caliper

PROCEDURE

- Set up the balance on a flat, stable surface.
- Use the vernier caliper to measure the length  $l$  of one of the cylinders. Make the measurement several times, average the values and record the average.
- Repeat the length measurement for the other cylinders.
- Use the vernier caliper to measure the diameter  $D$  of the one of the cylinders in several places. Average the values and record the average.
- Repeat the diameter measurement for the other cylinders.
- Weigh each of the cylinders and record the weights,  $W$ .
- Calculate the radius  $r$  for each cylinder ( $r = 1/2 \cdot D$ )
- Calculate the volume  $V$  of each cylinder using the formula  $V = \pi \cdot r^2 \cdot l$
- Calculate the density  $d$  of each cylinder using the formula  $d = W/V$ .

