

Zeeman Effect Apparatus

- **Complete Equipment Set for Studying Light Emission in a Magnetic Field**
- **Includes Video Capture Equipment and Computer Analysis Program**
- **Magnet Allows Transverse and Longitudinal Viewing**

The Zeeman Effect Apparatus examines the effect of a strong magnetic field on the green emission line of mercury at 545.1 nm wavelength. A magnetic field splits the degenerate 7s and 6p levels into three and five levels respectively, giving nine allowed transitions. The energy shifts of the splitting are tiny compared to the transition energy, so a high resolution spectrometer is necessary to observe them at the magnetic field strengths achievable in the laboratory. This resolution is usually achieved by using a high quality Fabry-Pérot étalon. This combination of these requirements has often discouraged the provision of Zeeman effect apparatus in the student laboratory. This apparatus addresses this issue.


The emitted light is polarized, and the nature and direction of the polarization depends on the orientation of the propagation direction to the magnetic field. For light emitted transversely to the field, the polarization is linear (π), with the polarization direction parallel to the electric field for one group of three transitions and perpendicular to the field for the remaining two groups. The components of these groups are usually not resolved sufficiently for easy measurement in practice due to overlap of their interference fringes. Light emitted parallel to the magnetic field is circularly polarized (σ), with only two groups of transitions emitted, showing opposing directions of polarization.

The electromagnet is mounted on a swivel on top of a power supply that also powers the slim low pressure mercury discharge lamp fitted between the magnet's pole pieces. A removable steel rod fills a long hole through one pole for transverse viewing; for parallel viewing, the magnet and lamp are swiveled 90° and the rod is replaced by a quarter wave plate for observing the circular polarization.


The optical system contains four elements and an optical bench. A polarizer and a 40mm diameter lens are mounted in a common housing. The polarizer can be rotated by 360° and clamped in any orientation. It serves as an analyzer for detecting the polarization state of the emitted light. The lens focuses the light for the spectrometer. An interference filter selects the emission line at 545.1 nm, and a Fabry-Pérot étalon with a fixed spacing of 2.0 mm acts as the resolving element. The étalon can be rotated in its housing and has three fine adjustment screws for setting exact parallelism of the mirrors. The étalon mount is a heavy steel block attached to a base plate with leveling screws. The steel lathe-bed style optical bench is equipped with four riders and the support rods for the optical elements carry height setting collars for easy fine alignment.

A CCD video camera with a mounting rod, power supply and PCI video capture card are supplied for direct live observation of the interference fringes. This makes measurement of fringe shifts on the computer monitor much easier than by visual observation using a conventional telescope arrangement.


Video capture and image analysis software is included in the equipment set. In addition to providing a live image of the Fabry-Pérot fringe pattern for viewing the fringe splitting as the magnetic field is increased and verifying the polarization states of the lines, an image can be captured and stored for analysis. The analysis module allows circles to be drawn on the image to exactly coincide with the positions of the fringes by manually inserting three reference points for each circle. To reduce image noise and improve positioning accuracy, a grayscale filter can be applied to the image before analysis. The software records the positions and radii of the circles drawn and, after circles have been added for all of the components of three successive fringes, calculates the wavelength shift for correlation with the magnetic field strength.



The étalon generates an interference pattern of concentric rings which is observed by a video camera.



Interference patterns are observed when the Zeeman apparatus is connected to a computer with a video capture card.



The initial appearance of the pattern with no imposed magnetic field.



This pattern shows the full range of Zeeman-split lines emitted in the transverse direction.



Triple fringes are strong π -polarized components that are used for measurement.



Fringes arising from the σ -polarized components are noticeably less intense.

Specifications

Lamp and Magnet:

Low pressure mercury discharge lamp, approx. 6mm diameter x 60mm long
 Electromagnet with power supply in base, swivels 90° for transverse or longitudinal viewing
 Current: 0, 0.3–2A
 Field strength: 1 T (max)
 Includes removable steel rod in one pole piece and plug-in quarter wave plate

Dimensions: 14cm x 29cm x 35cm

Power: 110 VAC

Weight: 43 kg.

Optical system:

Diameter of optical elements: 40mm

Interference filter: Central wavelength: 546.1 nm
 Transmission bandwidth: < 10nm

Fabry-Perot etalon:

Quartz mirror plates with 2.0mm spacing
 Central wavelength: 589.3nm
 Resolution ($\lambda/\delta\lambda$): > 2 x 10⁶
 High reflection bandwidth: 100nm

Computer and software requirements:

Requires a Windows-based computer with one available PCI slot running Windows 7 or 8 with at least 64MB of RAM and up to 1GB of hard disk space (for storing images)

Item No.

Description

ZEA001

Zeeman Effect Apparatus